

RELACS — a modular software platform for closed-loop experiments



Jan Benda

Department Biologie II
Ludwig-Maximilian
Universität München





“Traditional” Experiments

1. A set of stimuli and a more or less fixed experimental protocol are prepared



“Traditional” Experiments

1. A set of stimuli and a more or less fixed experimental protocol are prepared
2. The recordings are done on a few cells
($>$ one week of experimental work)
 - the raw voltage trace is the only feedback
 - changes during the running experiment often involve manual manipulations
 \Rightarrow precious recording time is wasted



“Traditional” Experiments

1. A set of stimuli and a more or less fixed experimental protocol are prepared
2. The recordings are done on a few cells
($>$ one week of experimental work)
 - the raw voltage trace is the only feedback
 - changes during the running experiment often involve manual manipulations
 \Rightarrow precious recording time is wasted
3. The data are analyzed offline



“Traditional” Experiments

1. A set of stimuli and a more or less fixed experimental protocol are prepared
2. The recordings are done on a few cells
($>$ one week of experimental work)
 - the raw voltage trace is the only feedback
 - changes during the running experiment often involve manual manipulations
 \Rightarrow precious recording time is wasted
3. The data are analyzed offline
4. The stimuli and the protocol are modified



“Traditional” Experiments

1. A set of stimuli and a more or less fixed experimental protocol are prepared
2. The recordings are done on a few cells
($>$ one week of experimental work)
 - the raw voltage trace is the only feedback
 - changes during the running experiment often involve manual manipulations
 \Rightarrow precious recording time is wasted
3. The data are analyzed offline
4. The stimuli and the protocol are modified
5. A new set of recordings is made





Closed-Loop Experiments

1. Present a stimulus
2. Record the response
3. Immediately analyze and visualize the data
4. Generate the next stimulus

Input:
Stimulus



Output:
Neuronal response



Output:
Next stimulus




Input:
Neuronal response



Closed-Loop Experiments

RELACS is designed as a framework for closed-loop experiments that

- considerably speed up this traditional approach
- offer novel experimental possibilities





Simple Closed-Loop Experiments

- Online visualization of processed data:
 - General infos, e.g. quality of spike detection, sensitivity of the cell, temperature, condition of animal, ...
 - Specific results, e.g. spike raster, firing rates, spike-triggered averages, ...
- ⇒ Speeds up manual closed-loop



Simple Closed-Loop Experiments

- Online visualization of processed data:
 - General infos, e.g. quality of spike detection, sensitivity of the cell, temperature, condition of animal, ...
 - Specific results, e.g. spike raster, firing rates, spike-triggered averages, ...
- ⇒ Speeds up manual closed-loop
- Set stimuli relative to the neuron's dynamic range
- Automatically control motorized electrodes
- ...



Advanced Closed-Loop Experiments


New experimental designs are possible:

- Optimal search for a neuron's receptive field.
- Search for stimuli that drive a neuron in an "optimal" way.
- Find set's of stimulus parameter that result in the same response (iso-response method).
- ...

Benda et al. (2007): "From response to stimulus: adaptive sampling in sensory physiology." *Curr. Opin. Neurobiol.* **17**: 430–436.



Example: Optimal Stimulus Ensembles



Machens et al. (2005) *Neuron* 17: 47–56.




Closed-Loop Experiments

In *RELACS* the closed-loop cycle can be freely programmed as a C++ plugin (*“experimental protocol”*).

The experimental-protocol plugins

- take recorded and pre-analyzed data
- perform analysis & display results
- generate next stimulus






Modular Design


RELACS core with flexible C++ Plugins for

- hardware abstraction
- data pre-processing (filter, spike detectors)
- passive and active control
- model





Screen Shot






Hardware Independent Protocols

RELACS integrates all hardware components.

Experimental protocols for *RELACS*

- are implemented independently of specific hardware
- can be used on all the different experimental setups in your lab without any modifications
- can be **shared** with other labs





Experimental Protocol Example

```
int Example::main( void ) {  
    // some initialization ...  
    OutData signal;  
    signal.setTrace( "LeftSpeaker" );  
    signal.sineWave( frequency, duration, amplitude );  
    SampleDataD rate( 0.0, duration, 0.001 );  
    for ( int counter=0; counter<Repeats; counter++ ) {  
        write( signal );  
        sleep( duration + pause );  
        EventData spikes( events( SpikeEvents[0] ),  
                           events( SpikeEvents[0] ).signalTime(),  
                           events( SpikeEvents[0] ).signalTime() + duration );  
        double meanrate = spikes.rate( 0.3*duration, duration );  
        spikes.addRate( rate, counter, GaussKernel( sigma ) );  
        P.clear();  
        P.plot( rate, 1000.0, Plot::Yellow, 2, Plot::Solid );  
        P.draw();  
        if ( meanrate < targetrate ) {  
            amplitude *= 2.0;  
            signal.sineWave( frequency, duration, amplitude );  
        }  
    }  
}
```



Why C++

- well structured (object oriented)
- platform independent
- efficient and controllable memory usage
- very fast



C++ Library for Data Analysis

Data structures (classes, container):

- *Array* — Basic 1-D vector
- *SampleData* — 1-D data vector with regularly sampled time axis
- *Map* — Sequence of $x|y$ data pairs

Algorithms:

- basic statistics (moments, quartiles, histogram)
- power spectra, coherence, transfer function
- linear fits
- non-linear fits (Simplex, Levenberg-Marquardt)



C++ Library for Data Analysis

Data structures (classes, container):

- *EventData* — Spikes and other point process data
- *EventList* — Multi-trial spike trains

Algorithms:

- firing rates (mean, PSTH binned/kernel, $1/ISI$)
- CV, Fano factor, ISI correlation
- vector strength, reliability, jitter
- mutual information (lower and upper bound)



Free and Open Source Software

RELACS is open source and free software distributed under the GNU General Public License (GPL).

- No hassle with licenses of commercial software.
- Add whatever new feature you need directly to the program.
- Share the program and your specific experimental protocols with your collaborators.
- Know what the data-analysis algorithms are doing!



Meta-data


Meta-data (“data about data”) describe the context of the raw data.

Meta-data are important for:

- Data management
- Data retrieving
- Data analysis
- Data sharing



Meta-Data Acquisition by RELACS





Meta-Data Acquisition by RELACS

RELACS - Relaxed Electrophysiology

File Plugins Devices Detectors/

Spikes-1

Minimum threshold 14.5
Threshold 36
Ratio threshold / size 40
Spike size 84.7

mV 0 50 100

Info Spectrum

Threshold [dB] 100 80 60 40 20
Frequency [kHz] 0 10 20 30 40
Best left 6.0

Firing rate [Hz] 10 5 0
Intensity [dB SPL] 50 100
Best 45.6

Search (ESC)

Stop Session Dialog

Meta Data Setup

Recording

Recording quality Good
Comment
Scientist Jan Benda
Temperature 22 °C
Humidity 54 %

Cell

Cell type Low-frequency receptor
Recording location Auditory nerve
Side Left
Depth 0 µm

Subject

Species Locusta migratoria
Sex Female
Age Middle-aged
Preparation in vivo dorsal

Save Discard Reset Cancel

on 0.9,6 (11/24/C)

F-I Field Search

Intensity [dB SPL] 40 60 80 100

Help (F3)

00:01 2008-12-02-ab/ 12% Frequency 4 kHz, Increment 4, Intensity 55 dB SPL, Loop 2



Meta-Data Acquisition by RELACS

RELACS records many meta-data:

- General infos about the experiment
(from the dialog)
- Main characteristics of the recorded cell
- All *RELACS*-controlled hardware settings
(e.g. sampling rate)
- All settings and version numbers of the
experimental protocols



The Meta-Data Problem

Name-value (+unit) pairs for:

- Stimuli
- Experimental settings
- Cell, preparation, experimental subject
- Hardware properties
- Analysis parameter
- etc.



The Meta-Data Problem

Name-value (+unit) pairs for:

- Stimuli
- Experimental settings
- Cell, preparation, experimental subject
- Hardware properties
- Analysis parameter
- etc.

But:

- What name to choose?
- What does it mean?
- How to share meta-data?



A Proposal

We need some kind of standard!

- flexible and extensible
- downloadable and parseable for seamless integration in software



A Proposal

We need some kind of standard!

- flexible and extensible
- downloadable and parseable for seamless integration in software

An XML schema explaining

- definitions
- dependencies
- meaning

of names (“tags”).



A Proposal

We need some kind of standard!

- flexible and extensible
- downloadable and parseable for seamless integration in software

An XML schema explaining

- definitions
- dependencies
- meaning

of names (“tags”).



In collaboration with G-node and LabLog by **Jan Grewe**
(lablog.sourceforge.net)

relacs_____

... enjoy your recordings



- Closed-loop experiments
- Hardware independent
- Data analysis libraries
 - Meta-data storage
 - Simulation mode
 - Dynamic clamp
- Open source, GPL, Linux
- ~ 120 000 lines of C++ code

by Jan Benda

[**www.relacs.net**](http://www.relacs.net)